

## How to add tutorial videos into ProSize data analysis software

To add the tutorial videos into ProSize data analysis software please follow these steps:

1. On the PC that has ProSize data analysis software installed go to (C:) drive
2. Click on the ProSize data analysis software folder
3. Right click and create a subfolder, name it as "Video"
4. Add/drag the downloaded "How-to videos" (.mp4 files, unzipped) into this folder.
5. Open ProSize data analysis software



6. Click on the "Help" menu
7. The How-to videos will be available under "Tutorial Videos"

