


# The Agilent ADM Flow Meter

Make productivity happen with this revolutionary flow meter—recalibrate right in your lab



## Key Features of the Agilent ADM Flow Meter




**Can record up to four flows** in volumetric and mass flow mode



**Is gas-type independent:** measures volumetric gas flow, no need to know the gas type composition



**Enables multi-tasking:** large OLEDdisplay lets you see the flow parameters without having to stand next to the device



**Instant replacement** of the new calibrated cartridge






**Direct download of new capabilities:** no need to ship flow meter for firmware upgrades



**Alert system:** warns you when the calibration cartridge needs replacing, minimizing the risk of volumetric measurement drift

## How the ADM Flow Meter Can Significantly Lower Your Ownership Costs

Annual Calibration Considerations	Traditional Flow Meter	Agilent ADM Flow Meter
 Shipping/paperwork charges	Yes	No
 Downtime	Yes	No
 Risk of damage during shipping	Yes	No



up to 45% saving on annual calibration costs\*

\*Based on average repair cost vs new cartridge

### Ease of use



One single Reference to order replacement cartridge



Easy to install/self installation of the cartridge

### Less Downtime



Instant replacement of the new calibrated cartridge

















No need for a back-up device



Replacement cartridge can be stored up to 1 year before installation



No failed recalibration

Annual Recalibration Costs	Conventional Flow Meters	ADM Flow Meter
 Cartridge recalibration	 \$165	 \$200
 Cartridge shipping <sup>1</sup>	 \$25-\$75*	\$0
 Repairs <sup>1</sup>	 \$125	\$0
 Backup meter during recalibration <sup>1</sup>	 \$73	\$0
 Total (Including recalibrated cartridge)	 \$388-\$438	 \$200
 Total added fees over 5 years	 \$1,115-1,365	\$0

\*Overseas shipment may be higher

<sup>1</sup>Annual added fees

See how the Agilent ADM Flow Meter makes it easier than ever to maintain NIST-Traceable standards at [www.agilent.com/chem/admflowmeter](http://www.agilent.com/chem/admflowmeter)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2018  
Published in the USA, September 1, 2018  
5994-0265EN